

Multiplace Hyperbaric Chamber Delivery Systems

Dick Clarke, CHT

The Multiplace Hyperbaric Chamber

Configuration & Operational Overview

Primary Training in Hyperbaric Medicine

Columbia, South Carolina

Operational, Research & Clinical Settings

- Military (diving & aviation)
- Commercial, professional & industrial diving worksites
- Recreational diving destination support
- Civil engineering (caisson & tunnel projects)
- Academic hyper-hypobaric programs
- Clinical hyperbaric medicine

The Class A Hyperbaric Chamber

'Human, multiple occupancy' NFPA-99 classification system

Multiplace hyperbaric delivery system fundamentals

- Chamber configuration
- Air compression & medical gas delivery systems
- Fire suppression capabilities
- Standard operating, clinical, safety & compliance characteristics

Chamber Nomenclature

- Main compartment vs. main lock/inner lock
- Entrance compartment vs. entry/outer lock
- Hatch vs. man way; 'door' on modern clinical units
- Medical lock vs. supply lock
- View port vs. window
- Inside attendant vs. tender/medic






Hatches open into compartment they seal


Outer compartment cannot be pressurized to greater degree than inner compartment

Hatches sealed by pressure differential, not by securing devices



The image shows two views of a hyperbaric chamber. On the left, the interior of the chamber is visible, showing a patient lying on a gurney. On the right, the exterior of the chamber is shown with the hatch open, revealing the interior. The exterior is white and has the letters 'HPO' and '6+2' printed on it, along with the 'ANST' logo.

Anti-suction protection, isolation valves; identification



The image shows a close-up of the anti-suction protection and isolation valves on a hyperbaric chamber. The valves are labeled with identification numbers and are connected to various tubes and hoses. The background is a white wall with some electrical wiring and components.

Chamber Air Compression Options

Low pressure compressor > LP volume receiver
125-200 psig

High pressure compressor > HP storage receiver
>2,000 psig

High pressure storage cylinders
>2,000 psig

Low Pressure Compressor Types

Internal combustion engine

Rotary screw



Compressor supplies chamber via air 'receiver'

small reserve for initial compression

receiver air initially cooler than direct compressor air



Air Purity

Compressors: responsibility of end users

Cylinders: responsibility of product manufacturer/supplier

~wide range of standards

Oxygen percent

Carbon dioxide; carbon monoxide

Oil vapor; water vapor

Sulphur dioxide; nitrogen dioxide

Nitric oxide, other potential contaminants

Filtration of Compressed Air Supply

Filter types

Particulate; to remove particles & dust

Activated carbon; removes odors & gases

Coalescing; aerosols to droplets, for drainage



Medical Gases

100% oxygen, clinical HBO therapy

Nitrox (nitrogen oxygen mixtures)

Heliox (helium oxygen mixtures)



Patient Oxygen Delivery Systems

- No ¼ turn valves on HP system; metering types only
- Pressure reduction to source, if HP cylinders
- Control panel O₂ analyzer
- 'Overboard dump' mandatory




Fire Suppression/Extinguishment Options

- Fire blankets
- Handheld extinguisher; water or foaming agents
- Water supplied hand line (hose)
- Water supplied deluge system

Fire Blankets

- Small portable deck decompression chambers
- Not permitted in Class A chambers (NFPA 99)



Handheld Fire Extinguisher




Mandatory Class A chambers (NFPA 99)

- Water & power (air pressure) must be independent of water deluge system
- Water pressure must be 50 psig > chamber working pressure
- One hose in every compartment, two in larger compartments
- Sufficient water volume for not < than 4 min



Water Deluge System

Mandatory Class A chambers (NFPA 99)



Water storage tank



Air pressure volume tank

Water Deluge System

Water & power (air pressure) supplies must be independent of headline

Water deluge must be delivered three seconds of activation



should be sufficient available volume to flow for one minute

should have sufficient stored power (pressure) to operate for 15 seconds

Water Deluge System

Manual activation & deactivation at control panel & inside chamber


Automatic activation not required; surveillance detection is (CCTV)


Evolution of flame detection - automatic activation standard

First edition NFPA-56 DT 1968

"A fixed fire extinguishing system shall be installed in all Class A chambers. It will be capable of manual, as well as automatic, activation"



NFPA 99, Chapter 19 1993

While the requirement for a fixed extinguishing system remains...

"Automatic fire detection systems are optional"

If installed...

"Surveillance fire detectors responsive to the radiation from flame shall be employed"

NFPA 99, Chapter 20, 2002

"Automatic fire detection systems shall not be required"

Inside Attendant

"Fitness to work in pressurized settings" along lines of HBO pts. but compressed air issue

Pregnancy temporary disqualification screening standards not well established

Nitrogen narcosis non-issue with routine HBO therapy complicates DCI cases tx > 100 fsw

historic standard of 165 fsw for CAGE & DCS

issues of personal safety, pt care & clinical decision-making

UHM JOURNAL, Vol. 2 - Policy Guidelines for Hyperbaric (October 2018) 673-674

Policy

UHM's Medical Fitness to Work Guidelines for Hyperbaric Inside Attendants

This document was written and reviewed by members of the Safety Committee (SA 10000) and approved by the UHM's Board of Directors

Purpose

The publication does not replace or supersede any other applicable regulations, standards, or policies. It is a guideline for reference only. It is not intended to be used as a legal document. The publication is subject to change without notice. The publication is the property of UHM and is not to be distributed outside of UHM.

The Board of Directors of the University of Maryland System (UMS) has approved this policy. It is a guideline for reference only. It is not intended to be used as a legal document. The publication is the property of UHM and is not to be distributed outside of UHM.

Authorship

This document is owned by the University of Maryland System. All rights reserved. No part of this document may be reproduced without the prior written permission of the University of Maryland System.

Alleman T. UHM 2018;45(2):231-247

Same pressure (barotrauma) related risks as patients

Also, at risk for decompression sickness several cases annually; two nurse fatalities

Procedures to limit DCS incidence


well hydrated; avoid cramped positions

O2 breathing interval prior to & during ascent

rotating personnel

avoid repetitive exposures

interval before flying >24 hrs.



UHM 2017, Vol. 44, No. 6 - HEALTH CARE WORKERS 673-607

Health care worker decompression sickness: incidence, risk and mitigation

Background: DCS is a rare but potentially fatal condition. Health care workers (HCWs) are at risk for DCS when they are exposed to high pressures during hyperbaric oxygen (HBO) therapy. This article reviews the incidence, risk factors, and mitigation strategies for DCS in HCWs.

Abstract

Decompression sickness (DCS) is a rare but potentially fatal condition that can affect anyone who is exposed to high pressures during hyperbaric oxygen (HBO) therapy. Health care workers (HCWs) are at risk for DCS when they are exposed to high pressures during HBO therapy. This article reviews the incidence, risk factors, and mitigation strategies for DCS in HCWs.

Keywords

decompression sickness; hyperbaric oxygen; health care workers; incidence; risk; mitigation

Clarke R. UHM 2017;44(6):509-519

Routine Chamber Operations

Chamber compression on air to pre-determined depth
*traditionally 45 fsw/2.36 ATA; 33 fsw/2.0 ATA common with hoods
mandatory BIBS for every occupant
O₂ breathing upon arrival tx pressure*

Chamber O₂ constantly monitored; upper limit 23.5%
*most common source of O₂ leaks is BIBS
air flushing to maintain acceptable range*

Entry lock at 1.0 ATA

Medical lock for small supplies transfer
ampules; vials



Patient/Chamber Safety Perspectives

No direct patient grounding required (as required for monoplace)
exception when O₂ atmospheres > 23.5% employed

Chamber grounding per monoplace, plus internal conductive surfaces
enclosing electrical circuits

Ear protection during compression & ventilation

Chamber pressure relief valves
*one for chamber max working pressure
optional second for max O₂ breathing pressure*

Battery powered equipment
*sealed/pressure resistant; no in-chamber charging; no battery changing
personal items (cell phone; laptop; tablet; pager; entertainment) prohibited*

"Intrinsically Safe" Hyperbaric Practice

Keeping level of electrical energy too low to cause ignition
thereby preventing sparks & keeping temperature low

Device designs that exclude oxygen
plus, purging device with inert gas

Device strong enough to contain explosion

Moving devices outside hazardous (chamber) area

